

FLANGE SIZE


can damage your milk supply + breast tissue



flange too big too much areola pulled inside tunnel, nipple base not touching walls



flange perfect fit nipple fits comfortably inside, glides along walls



flange too small only tip of nipple pulled inside tunnel, squeezed tight against walls



Measure nipple diameter only, focusing on the base.





Breast shield (flange) should measure 2-4mm larger than nipple measurement. This can vary.



nipple size:













24mm





recommended starting flange size: 13mm 15mm

17mm

19mm

21mm

27mm

30mm

1 . 1